| AUII/D. | Environmental Analysis Teaching | Date: 1/27/2024 | SOP No. 37b v0.92 |
|-------------------|---------------------------------|---------------------------|-------------------|
| POMONA | and Research Laboratory | | |
| | Standard Operating Procedure | Title: Compost Maturity | |
| POMONA COLLEGE | Approved By: Los Huertos | Revised: January 28, 2024 | |

1. Scope and Application

- 1.1 The Solvita Compost Maturity Test is a unique procedure using two test probes to simultaneously measure CO2 and NH3, the two most prominent gases indicating activity and stability of compost.
- 1.2 This test is widely-recognized worldwide for validating compost maturity.
- 1.3 This easy-to-use application requires a small sample placed in the incubation jars and the results are read after 4-hrs of exposure. The basic kit offers a Solvita color chart to determine the level of activity from which a Maturity Index is calculated using guidelines provided in the manual.
- 1.4
- 1.5 Numerous studies... (Vargas et al., 2005).

2. Acknowledgements

2.1 This SOP is based on the Solvita Field Test Kit Manual (Solvita, 2014), version 8.

3. Definitions

Compost Compost Maturity Maturity Index

4. Biases and Interferences

- **4.1** Biases and interferences can come from...
- 5. Health and Safety
- 6. Personnel & Training Responsibilities

Students using this SOP should be trained for the following SOPsa:

- SOP1
- SOP2

Author: Marc Los Huertos |
File: SolvitaCompostMaturity.tex

Page: 1 of 5

SOP: 37b v0.92 (Revised: January 28, 2024)

7. Required Materials

- Test Soils
- # XX Seives (3/8" or 10 mm)
- Solvita Jars
- Individually wrapped CO2 & NH3 Probes (must remain refrigerated)
- DCR Field Test Unit (with 9 volt battery)
- AWS digital scale (field) or Ohaus digital scale (lab)

8. Estimated Time

8.1 This will take 20 minutes to collect and prepare compost sample. An equalibrium stage ranges from 1 to 24 hours depending on conditions. The incubation requires 4 hours. Finally, an additional 10 minutes is needed to read the probes for a total of 5.5 to 28.5 hours.

9. Sample Collection, Preservation, and Storage

Obtain and Prepare Sample

- 9.1 Take several grab samples to prepare a compost by mixing all sub-samples resprese-native of the enre compose. Remove larg weed chips and other objects. A 3/8" (10 mm) seive is recommended before laoding jar.
- **9.2** Check moisture content of sample using squeeze test. Squeeze a handful of compost and release. Water should appear between fingers, but not drip out. Adjust water content if needed. Allow compost to equilibrate for 24 hours.
- **9.3** If compost is warm or frozen, let it equilibrate to room temperature for 24 hours before testing.

10. Procedure

- 10.1 Seive compost to remove large particles.
- 10.2 Carefully fill the Solvita jar to the fill lines. To get the proper density tap the jar gently while filling to fill line. Optionally, the proper weight in grams per jar corresponding to actual field density is found in the Table 1.
- 10.3 Let the sample "air-out" briefly in the jar for up to one-hour prior to starting the test. If the sampl was take directly from a very hot or frozen pile, it is advisable to allow it to stand for 24 hours before starting the test.
- 10.4 The Solvita maturity method requires two tests carried out simultanious in the same 4-hour period. Terr open both the pouches marked "High CO2" and "Ammonia"

Author: Marc Los Huertos Page: 2 of 5

SOP: 37b v0.92 (Revised: January 28, 2024)

Table 1: Compost weight per jar adjusted for actual field density.

| lbs/yr3 | kg/m3 | g/jar |
|---------|-------|-------|
| 500 | 300 | 30 |
| 800 | 475 | 50 |
| 1,000 | 600 | 60 |
| 1,200 | 700 | 70 |

and carefully remove each probe and avoid touching the gel surface. The gel in the proble is color-coded: the carbon dioxide probe is purple and the ammonia probe is yellow. Once the pouch is open, the test should be started immediately.

- 10.5 Carefully insert both probes into the compost without allowing compost to touch the gel surface. Push the probes so the bottom of the probe hits the jar's bottom. Do not jostle the composite or allow dust to coat the surface of the probe.
- 10.6 Screw on lid tightly and allow to incubate out of direct sunlight and wait 4 hours. Record temperature and try to keep the jars at a constant temp for the duration of the test, between 20-25°C (68-77°F).
- 10.7 Remove lid after 4 hours and remove the probes.
- 10.8 Read the probes immediately after removing from the jar.

Reading the Results

- 10.9 Turn on DCR Field test unit and insert probe to get CO2 and NH3 color. The probe must go into the DCR with gel side up press the read button.
- 10.10 Insert the NH3 probe into the DCR and press the read button. Compare color to the visual color key. Record the NH3 reading.
- 10.11 If the DCR does not work, Compare the color of the gel to the color chart, finding the closest match. Record the CO2 and NH3 readings.

Clean Up

- **10.12** Discard compost sample and probe in an appropriate manner.
- 10.13 Wash the Solvita jar with water and dry for reuse. Label jar top with event to track the number of times the jar has been used. Disscard the jar after 4 uses.

11. Data Analysis and Calculations

11.1 Calculate the maturity index using the following formula:

Author: Marc Los Huertos

SOP: 37b v0.92 (Revised: January 28, $\overline{2024}$)

$$MI = \frac{CO2 \times 0.5 + NH3 \times 0.5}{100} \tag{1}$$

Compost Maturity Index

11.2 See Figure! 1 below for interpretation:

12. Applications of Results

Status and Conditions of Compost Processes

Managing Aeration Sufficiency

Selecting Best Use of Composed Based on Maturity

Ammoiniua Emissions of Mantures & Composts

13. QC/QA Criteria

- 13.1 The probe show Lot No. and Expiration Dates on the package. Solvita kits are pre-calibrated and packaged for the highest quality prior to shipping. The sealed probes should be the "Control Color" when the foil pack is opened (see color chart). Shelf-life is improved by refrigeration. Do not allow gels to freeze.
- **13.2** The plastic jars may be reused 4 times, the discarded.
- **13.3** If the foil packs are damaged or the jar is cracked then the test may not work properly.
- **13.4** A sample that is too dry may give a false postive maturiety test (although goty ammonia is still volatile).

14. Trouble Shooting

Table 2: Trouble Shooting Guide

| Problem | Cause | Solution |
|-------------------------------|--------------------|--------------------------|
| Compost is fresh but test re- | Compost is too dry | Add water to compost and |
| sults indcate "mature" | | retest |

Author: Marc Los Huertos Page: 4 of 5

SOP: 37b v0.92 (Revised: January 28, 2024)

15. References

References

Solvita, L. (2014). Official Solvita Guideline: Compost Emissions Test: A New Tool for Compost Maturity Assessment.

Vargas, D. C., Monedero, M. S., Urpilainen, S., Kamilaki, A., and Stentiford, E. (2005). Assessing the stability and maturity of compost at large-scale plants. *Ingeniería*, 9(2):25–30.

Author: Marc Los Huertos Page: 5 of 5