

	Environmental Analysis Teaching and Research Laboratory	Date: 1/27/2024	SOP No. 37b v0.92
	Standard Operating Procedure	Title: Compost Maturity	
	Approved By: Los Huertos	Revised: January 28, 2024	

1. Scope and Application

1.1 The Solvita Compost Maturity Test is a unique procedure using two test probes to simultaneously measure CO₂ and NH₃, the two most prominent gases indicating activity and stability of compost.

1.2 This test is widely-recognized worldwide for validating compost maturity.

1.3 This easy-to-use application requires a small sample placed in the incubation jars and the results are read after 4-hrs of exposure. The basic kit offers a Solvita color chart to determine the level of activity from which a Maturity Index is calculated using guidelines provided in the manual.

1.4

1.5 Numerous studies... ([Vargas et al., 2005](#)).

2. Acknowledgements

2.1 This SOP is based on the Solvita Field Test Kit Manual ([Solvita, 2014](#)), version 8.

3. Definitions

3.1 Compost

3.2 Compost Maturity

3.3 Maturity Index

4. Biases and Interferences

4.1 Biases and interferences can come from...

5. Health and Safety

6. Personnel & Training Responsibilities

Students using this SOP should be trained for the following SOPsa:

- SOP1
- SOP2

7. Required Materials

- Test Soils
- # XX Sieves (3/8" or 10 mm)
- Solvita Jars
- Individually wrapped CO₂ & NH₃ Probes (must remain refrigerated)
- AWS digital scale (field) or Ohaus digital scale (lab)

8. Estimated Time

8.1 This will take 15 minutes in the lab and 24 hours for the Solvita Jars to incubate.

9. Sample Collection, Preservation, and Storage

Obtain and Prepare Sample

- 9.1 Take several grab samples to prepare a compost by mixing all sub-samples representative of the entire compost. Remove large weed chips and other objects. A 3/8" (10 mm) sieve is recommended before loading jar.
- 9.2 Check moisture content of sample using squeeze test. Squeeze a handful of compost and release. Water should appear between fingers, but not drip out. Adjust water content if needed. Allow compost to equilibrate for 24 hours.
- 9.3 If compost is warm or frozen, let it equilibrate to room temperature for 24 hours before testing.

10. Procedure

- 10.1 Place a clean Solvita Jar on the scale and tare the weight of the jar
- 10.2 Add 100 ± 5 grams of soil using the fill line as a guide.
- 10.3 Unwrap and place CO₂ probe in jar NOTE: Handle the probe only by the handle avoid anything from touching the gel surface.
- 10.4 Screw on lid tightly– and wait 24 hours. Record temperature and try to keep the jars at a constant temp for the duration of the test
- 10.5 Remove lid after 24 hours

Reading the Results

10.6 Turn on DCR Field test unit and insert probe to get CO₂ color. The probe must go into the DCR with gel side up press the read button. Compare color to the visual color key. Record the CO₂ reading.

10.7 Insert the NH₃ probe into the DCR and press the read button. Compare color to the visual color key. Record the NH₃ reading.

11. Data Analysis and Calculations

11.1 Calculate the maturity index using the following formula:

$$MI = \frac{CO_2 \times 0.5 + NH_3 \times 0.5}{100} \quad (1)$$

11.2 See Table 1 below for interpretation:

12. QC/QA Criteria

12.1 Solvita kits are pre-calibrated and packaged for the highest quality prior to shipping. The sealed probes should be the "Control Color" when the foil pack is opened (see color chart). If the foil packs are damaged or the jar is cracked then the test may not work properly. The probe show Lot No. and Expiration Dates on the package. The plastic jars may be reused 4 times, the discarded. Shelf-life is improved by refrigeration. Do not allow gels to freeze.

13. Trouble Shooting

14. References

References

Solvita, L. (2014). *Official Solvita Guideline: Compost Emissions Test: A New Tool for Compost Maturity Assessment*.

Vargas, D. C., Monedero, M. S., Urpilainen, S., Kamilaki, A., and Stentiford, E. (2005). Assessing the stability and maturity of compost at large-scale plants. *Ingeniería*, 9(2):25–30.

14.1 APHA, AWWA. WEF. (2012) Standard Methods for examination of water and wastewater. 22nd American Public Health Association (Eds.). Washington. 1360 pp. (2014).